

Product Datasheet

ACTA2 antibody (FITC) (orb466933)

Description	ACTA2 antibody (FITC)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	FC, IF
Immunogen	KLH conjugated synthetic peptide derived from human Actin alpha
Target	alpha smooth muscle Actin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P62736
Dilution Range	Flow-Cyt=1ug/Test,IF=1:500-1000
Expiration Date	12 months from date of receipt.