



Product Datasheet

CRKL antibody (orb46748)

Description	CRKL antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Crk-like protein (98-303AA)
Target	CRKL
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P46109
Dilution Range	WB:1:2000-10000, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2(lane 1), ...



Immunofluorescence analysis of HeLa cell...