


Product Datasheet

ETS2 antibody (orb46768)

Description	ETS2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Protein C-ets-2 protein (9-185AA)
Target	ETS2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P15036
Dilution Range	WB:1:500-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of K562(lane 1), H...