


Product Datasheet


SNCA antibody (orb46816)

Description	Goat polyclonal antibody to SNCA
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	Alpha-synuclein
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Concentration	0.5 mg/ml
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 17kDa band observed in Human Brain (Amygdala) and fetal Mouse Brain lysates (calculated MW of 14.5kDa according to Human NP_000336.1 and Mouse NP_033247.1). This molecular weight is routinely observed by other sources. Recommended concentration: 0.. IHC: Paraffin embedded Human Brain (Cerebellum). Recommended concentration: 5µg/ml.
Protein Sequence	ATGFVKKDQLGK
Clonality	Polyclonal
MW	14.5
Entrez	6622
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis
of Human Amygdala

...



Immunohistochemical
staining of Human
Ce...