

## Product Datasheet

### Thyroxine Binding Globulin antibody (HRP) (orb468792)

**Description**

Rabbit polyclonal antibody against

conjugated to HRP

**Species/Host**

Rabbit

**Conjugation**

HRP

**Tested**

ELISA

**Applications**
**Immunogen**

KLH conjugated synthetic peptide  
derived from human Thyroxine  
Binding Globulin

**Target**

SERPINA7

**Preservatives**

0.01M TBS(pH7.4) with 1% rAlbumin,  
0.03% Proclin300 and 50% Glycerol.

**Form/Appearance**

Lyophilized or Liquid

**Concentration**

- 100 µg (in 100 µl): 1 mg/ml

**Storage**

Store at -20 °C for one year. Avoid  
repeated freeze/thaw cycles. The  
lyophilized antibody is stable at room  
temperature for at least one month  
and for greater than a year when kept  
at -20°C. When reconstituted in sterile  
pH 7.4 0.01M PBS or diluent of  
antibody the antibody is stable for at  
least two weeks at 2-4 °C.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**MW**

44kDa

**Uniprot ID**
**P05543**
**Expiration Date**

12 months from date of receipt.