

Product Datasheet

IL-1 Alpha Propeptide antibody (HRP) (orb469950)

Description	IL-1 Alpha Propeptide antibody (HRP)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	HRP
Tested Applications	ELISA, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from mouse interleukin-1 alpha Propeptide
Target	IL-1 Alpha Propeptide
Preservatives	lyophilized from a solution containing protein, potassium phosphate buffer and sodium chloride. No preservatives are added.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P01583
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.