

Product Datasheet

PLEKHO1 antibody (Biotin) (orb47090)

Description	PLEKHO1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Pleckstrin homology domain-containing family O member 1 protein (143-404AA)
Target	PLEKHO1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q53GL0
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.