



Product Datasheet

Prostaglandin E Receptor EP2 antibody (orb4710)


Description	Prostaglandin E Receptor EP2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human Prostaglandin E Receptor EP2 (271-358/358 aa)
Target	PTGER2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40 kDa
Uniprot ID	P43116
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Spleen (Mouse) ...



Flow cytometric analysis of A549 cells U...



Western blot analysis of Placenta (Mouse...