



Product Datasheet

GRAMD4 antibody (orb47115)


Description	GRAMD4 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human GRAM domain-containing protein 4 protein (23-161AA)
Target	GRAMD4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6IC98
Dilution Range	WB:1:500-5000, IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Western blot aanalysis of Mouse brain ti...



Western blot aanalysis of MCF-7 whole ce...



Immunohistochemistry of paraffin-embedde...