

Product Datasheet

HBcAg antibody (HRP) (orb471195)

Description	HBcAg antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA, ICC
Immunogen	KLH conjugated synthetic peptide derived from Hepatitis B Virus Core Antigen
Target	HBcAg
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P03148
Dilution Range	ICC=1:100,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.