

## Product Datasheet

### KBTBD6 antibody (FITC) (orb47125)

|                        |                                                                                      |
|------------------------|--------------------------------------------------------------------------------------|
| <b>Description</b>     | KBTBD6 antibody (FITC)                                                               |
| <b>Species/Host</b>    | Rabbit                                                                               |
| <b>Reactivity</b>      | Human                                                                                |
| <b>Conjugation</b>     | FITC                                                                                 |
| <b>Immunogen</b>       | Recombinant Human Kelch repeat and BTB domain-containing protein 6 protein (1-139AA) |
| <b>Target</b>          | KBTBD6                                                                               |
| <b>Preservatives</b>   | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol                       |
| <b>Form/Appearance</b> | Liquid                                                                               |
| <b>Storage</b>         | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                        |
| <b>Note</b>            | For research use only                                                                |
| <b>Isotype</b>         | IgG                                                                                  |
| <b>Clonality</b>       | Polyclonal                                                                           |
| <b>Uniprot ID</b>      | <a href="#">Q86V97</a>                                                               |
| <b>Dilution Range</b>  | ELISA                                                                                |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                      |