

Product Datasheet

ANXA4 antibody (FITC) (orb47297)

Description	ANXA4 antibody (FITC)
Species/Host	Rabbit
Reactivity	Canine
Conjugation	FITC
Immunogen	Recombinant Dog Annexin A4 protein (2-319AA)
Target	ANXA4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P50994
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.