

Product Datasheet

Irf3 antibody (orb47427)

Description	Irf3 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Mouse Interferon regulatory factor 3 protein (1-419AA)
Target	Irf3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P70671
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.