

Product Datasheet

GFP-Tag antibody (HRP) (orb474838)

Description	GFP-Tag antibody (HRP)
Species/Host	Mouse
Conjugation	HRP
Tested Applications	ELISA, ICC, WB
Immunogen	Recombinant Green Fluorescence Protein
Target	GFP
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P42212
Dilution Range	WB=1:5000-20000,ICC=1:100-500,ELISA=1:5000-20000
Expiration Date	12 months from date of receipt.