

Product Datasheet

MLKL (phospho-Ser358) antibody (HRP) (orb474849)

Description

Mouse monoclonal antibody against

HRP

Species/Host

Mouse

Conjugation

HRP

Tested

ICC, IHC-Fr, IHC-P

Applications
Immunogen

KLH conjugated synthetic peptide derived from human phospho-MLKL (Ser358)

Target

MLKL

Preservatives

Constituents: 0.01M PBS, pH 7.4 with 10 mg/mL rAlbumin and 0.1% Gentamicin, 50% glycerol. Or Lyophilized. Buffer = 0.01M PBS, pH 7.4 with 10 mg/mL rAlbumin and 0.1% Gentamicin. Reconstitute with sterile distilled water.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Monoclonal

MW

54kDa

Uniprot ID

Q9D2Y4

Dilution Range

IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500

Expiration Date

12 months from date of receipt.