

Product Datasheet

CD34 antibody (HRP) (orb474860)

Description	CD34 antibody (HRP)
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	HRP
Tested Applications	ELISA, IHC-P, WB
Immunogen	Recombinant human CD34
Target	CD34
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P28906
Dilution Range	WB=,IHC-P=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.