

## Product Datasheet

### Cytokeratin 18 antibody (HRP) [Out of stock] (orb474998)

|                            |                                                                                                   |
|----------------------------|---------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Cytokeratin 18 antibody (HRP) [Out of stock]                                                      |
| <b>Species/Host</b>        | Rabbit                                                                                            |
| <b>Conjugation</b>         | HRP                                                                                               |
| <b>Tested Applications</b> | FC, ICC, IF, IHC-Fr, IHC-P, WB                                                                    |
| <b>Immunogen</b>           | Recombinant human Cytokeratin 18                                                                  |
| <b>Target</b>              | CK18                                                                                              |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.                             |
| <b>Form/Appearance</b>     | Liquid                                                                                            |
| <b>Concentration</b>       | 1mg/ml                                                                                            |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.                  |
| <b>Note</b>                | For research use only                                                                             |
| <b>Isotype</b>             | IgG                                                                                               |
| <b>Clonality</b>           | Monoclonal                                                                                        |
| <b>Uniprot ID</b>          | <b>P05783</b>                                                                                     |
| <b>Dilution Range</b>      | WB: 1:500-2000, IHC-P: 1:100-500, IHC-F: 1:100-500, ICC: 1:100, IF: 1:100-500, Flow-Cyt: 2ug/Test |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                   |