

Product Datasheet

Neuromedin U antibody (HRP) (orb475110)

Description

Rabbit polyclonal antibody against

Species/Host	Rabbit
Conjugation	HRP
Tested Applications	ICC, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from mouse Neuromedin U
Target	NMU
Preservatives	Constituents: 0.01M PBS, pH 7.4 with 10 mg/mL rAlbumin and 0.1% Gentamicin, 50% glycerol. Or Lyophilized. Buffer = 0.01M PBS, pH 7.4 with 10 mg/mL rAlbumin and 0.1% Gentamicin. Reconstitute with sterile distilled water.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20kDa
Uniprot ID	P48645
Dilution Range	IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.