

Product Datasheet

HP1 Gamma antibody (HRP) (orb475216)

Description	HP1 Gamma antibody (HRP)
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	HRP
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P, WB
Immunogen	Recombinant human HP1 gamma Protein
Target	HP1 gamma
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	21kDa
Uniprot ID	Q13185
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,IHC- F=1:100-500,ICC=1:100- 500,ELISA=1:1000-5000
Expiration Date	12 months from date of receipt.