

Product Datasheet

MCherry-Tag antibody (HRP) (orb475224)

Description	MCherry-Tag antibody (HRP)
Species/Host	Mouse
Reactivity	All
Conjugation	HRP
Tested Applications	WB
Immunogen	Recombinant Protein mCherry
Target	mCherry
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	25kDa
Uniprot ID	X5DSL3
Dilution Range	WB=1:1000-5000
Expiration Date	12 months from date of receipt.