

Product Datasheet

p53 (Acetyl K382) antibody (HRP) (orb475246)

Description	p53 (Acetyl K382) antibody (HRP)
Species/Host	Mouse
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised acetylpeptide derived from human P53 around the acetylation site of Acetyl K382
Target	p53 (Acetyl K382)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P04637
Dilution Range	IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.