

Product Datasheet

Beta-Actin antibody (HRP) (orb475250)

Description	Beta-Actin antibody (HRP)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	HRP
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P, WB
Immunogen	Synthetic MAP peptide derived from human beta-Actin
Target	beta-Actin (Loading Control)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P60709
Dilution Range	WB=1:5000-20000,IHC-P=1:100-500,IHC-F=1:200-1000,ICC=1:100,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.