

Product Datasheet

PRPF4 antibody (HRP) (orb475506)

Description

Rabbit polyclonal antibody against PRPF4

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	HRP
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from human PRPF4
Target	PRPF4B
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	58kDa
Uniprot ID	Q13523
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.