

Product Datasheet

ALPK3 antibody (HRP) (orb475662)

Description	ALPK3 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human ALPK3
Target	ALPK3
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96L96
Dilution Range	IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.