

Product Datasheet

CDX2 antibody (orb47567)

Description	CDX2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Homeobox protein CDX-2 protein (12-116AA)
Target	CDX2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q99626
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.