


Product Datasheet

CaMK2 alpha (phospho-Thr305) antibody (orb4757)

Description	CaMK2 alpha (phospho-Thr305) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human CaMK2 alpha around the phosphorylation site of Thr305 (IL(p-T)TM)
Target	phospho-CaMK2 alpha (Thr305)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	54 kDa
Uniprot ID	Q9UQM7
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot: Eye (Mouse); Lung (Mouse); Pr...