

Product Datasheet

TLR7 antibody (HRP) (orb47596)

Description	TLR7 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Toll-like receptor 7 protein (899-1022AA)
Target	TLR7
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NYK1
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.