

Product Datasheet

EGM_08965 antibody (Biotin) (orb47602)

Description	EGM_08965 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Macaca fascicularis Myelin basic protein (1-304AA)
Target	MBP
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	G7PWA0
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.