



Product Datasheet

LRP8 antibody (orb47635)

Description	LRP8 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Low-density lipoprotein receptor-related protein 8 protein (131-332AA)
Target	LRP8
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q14114
Dilution Range	WB:1:500-5000, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of SH-SY5Y(lane 1)...



Immunofluorescence analysis of HepG2 cel...