

Product Datasheet

ARHGAP11B antibody (FITC) (orb47645)

Description	ARHGAP11B antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Rho GTPase-activating protein 11B protein (99-231AA)
Target	ARHGAP11B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q3KRB8
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.