

Product Datasheet

Gfer antibody (Biotin) (orb47666)

Description	Gfer antibody (Biotin)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Rat FAD-linked sulfhydryl oxidase ALR protein (1-198AA)
Target	Gfer
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q63042
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.