



Product Datasheet

SLC16A9 antibody (orb47667)


Description	SLC16A9 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Monocarboxylate transporter 9 protein (185-304AA)
Target	SLC16A9
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7RTY1
Dilution Range	WB:1:500-5000, IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of SH-SY5Y(lane 1)...



Immunohistochemical staining of human sm...



Immunofluorescence analysis of PC-3 cell...