

Product Datasheet

ANKRD11 antibody (orb47671)

Description	ANKRD11 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Ankyrin repeat domain-containing protein 11 protein (1002-1276AA)
Target	ANKRD11
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6UB99
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.