

Product Datasheet

PPP1R13B antibody (FITC) (orb47717)

Description	PPP1R13B antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Apoptosis-stimulating of p53 protein 1 protein (375-675AA)
Target	PPP1R13B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96KQ4
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.