

Product Datasheet

Phosphotyrosine antibody (HRP) [Out of stock] (orb478236)

Description	Rabbit polyclonal antibody against Phosphotyrosine conjugated to HRP
Species/Host	Rabbit
Conjugation	HRP
Immunogen	KLH conjugated Phosphotyrosine
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
CAS Number	21820-51-9
Expiration Date	12 months from date of receipt.