

Product Datasheet

ABO antibody (orb47827)

Description	ABO antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Histo-blood group ABO system transferase protein (235-354AA)
Target	ABO
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P16442
Dilution Range	WB:1:200-2000
Expiration Date	12 months from date of receipt.