

Product Datasheet

Ghrelin 28 (n-octanoyl Ser26) antibody (HRP) (orb478287)

Description	Ghrelin 28 (n-octanoyl Ser26) antibody (HRP)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	HRP
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised n-octanoyl-peptide derived from human Ghrelin 28 around the n-octanoyl site of Ser26)
Target	Ghrelin 28 (n-octanoyl Ser26)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UBU3
Dilution Range	IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.