

Product Datasheet

AAV5 capsid protein VP1 antibody (HRP) (orb479018)

Description

Rabbit polyclonal antibody against

	HRP
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from Adeno-Associated Virus 5 capsid protein VP1 (46-80aa)
Target	AAV5
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	80kDa
Dilution Range	IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558