


Product Datasheet

c-Kit (phospho-Tyr703) antibody (orb4796)

Description	c-Kit (phospho-Tyr703) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human c-Kit around the phosphorylation site of Tyr703 (AL(p-Y)KN)
Target	Phospho-c-Kit (Tyr703)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	105 kDa
Uniprot ID	P10721
Dilution Range	Flow-Cyt=1µg /test
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
K562 cells
u...