



Product Datasheet

MST1R antibody (orb4801)


Description	MST1R antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human Ron (151-250/1400 aa)
Target	MST1R
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	30/119/150 kDa
Uniprot ID	Q04912
Dilution Range	Flow-Cyt=1ug/Test
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
A431 cells
u...



Flow
cytometric
analysis of
A431 cells
u...



Flow
cytometric
analysis of
A431 cells
u...