

Product Datasheet

Acetyl-Histone H3 (K23) antibody (HRP) (orb480131)

Description

Acetyl-Histone H3 (K23) antibody (HRP)

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

HRP

Tested

ELISA, ICC, IHC-Fr, IHC-P

Applications
Immunogen

KLH conjugated synthesised
acetylpeptide derived from human
Histone H3 around the acetylation site
of K23

Target

Histone H3 (Acetyl K23)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin,
0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Shipped at 4°C. Store at -20 °C for one
year. Avoid repeated freeze/thaw
cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

CAS Number

Histone H3.1

Uniprot ID
P68431
Dilution Range

IHC-P=1:100-500,IHC-F=1:100-
500,ICC=1:100-500,ELISA=1:5000-
10000

Expiration Date

12 months from date of receipt.