

Product Datasheet

EPO/Erythropoietin Antibody (orb48037)

Description

EPO/Erythropoietin Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

E. coli-derived human EPO recombinant protein (Position: M1-R166).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

21307 MW

Uniprot ID

P01588


Dilution Range

Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, By Heat ELISA , 0.1-0.5µg/ml, Human, -Western blot, 0.1-0.5µg/ml, Human


Expiration Date

12 months from date of receipt.

1 2 3



Western blot analysis of Lane 1: Recombi...



Immunohistochemical analysis of formalin...