

Product Datasheet

Interleukin-4 IL4 Antibody (orb48042)

Description

Interleukin-4 IL4 Antibody



Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

WB analysis
of IL-4 using
anti-IL-4
anti...

Immunogen

E. coli-derived human IL-4 recombinant protein
(Position: H25-S153).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a
concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After
reconstitution, at 4°C for one month. It can also be
aliquotted and stored frozen at -20°C for six months.
Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

By Heat: Boiling the paraffin sections in 10mM citrate
buffer, pH6.0, for 20mins is required for the staining of
formalin/paraffin sections. Other applications have not
been tested. Optimal dilutions should be determined by
end users. . Add 0.2ml of distilled water will yield a
concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

17492 MW

Uniprot ID

P05112

Dilution Range

Immunohistochemistry (Paraffin-embedded Section),
0.5-1µg/ml, Human, By Heat ELISA , 0.1-0.5µg/ml,
Human, -Western blot, 0.1-0.5µg/ml, Human

Expiration Date

12 months from date of receipt.