

Product Datasheet

Leptin antibody (HRP) (orb480466)

Description	Rabbit polyclonal antibody against Leptin conjugated to HRP
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human Leptin
Target	Leptin (HRP)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid: 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	16 kDa
Uniprot ID	P41159
Dilution Range	ELISA: 1:100-1000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.