



Product Datasheet

CEBPA antibody (orb48049)

Description	CEBPA antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human CCAAT/enhancer-binding protein alpha protein (1-124AA)
Target	CEBPA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P49715
Dilution Range	WB:1:200-3000, IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human pa...



Western blot analysis of A549(lane 1), H...