



Product Datasheet

NDEL1 (phospho-Ser242) antibody (orb4808)


Description	NDEL1 (phospho-Ser242) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human NDEL1 around the phosphorylation site of Ser242 (GT(p-S)PL)
Target	phospho-NDEL1 (Ser242)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	38 kDa
Uniprot ID	Q9GZM8
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ce...



Western blot analysis of mouse intestine...



Western blot analysis of mouse Epitheliu...