

Product Datasheet

NDEL1 (phospho-Thr219) antibody (orb4809)

Description	NDEL1 (phospho-Thr219) antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human NDEL1 around the phosphorylation site of Thr219 (PA(p-T)PV)
Target	phospho-NDEL1 (Thr219)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	38 kDa
Uniprot ID	Q9GZM8
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.