

Product Datasheet

EIF3J antibody (FITC) (orb48162)

Description	EIF3J antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Eukaryotic translation initiation factor 3 subunit J protein (5-115AA)
Target	EIF3J
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O75822
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.