

## Product Datasheet

### GNAQ antibody (FITC) (orb48254)

|                        |                                                                                           |
|------------------------|-------------------------------------------------------------------------------------------|
| <b>Description</b>     | GNAQ antibody (FITC)                                                                      |
| <b>Species/Host</b>    | Rabbit                                                                                    |
| <b>Reactivity</b>      | Human                                                                                     |
| <b>Conjugation</b>     | FITC                                                                                      |
| <b>Immunogen</b>       | Recombinant Human Guanine nucleotide-binding protein G(q) subunit alpha protein (1-359AA) |
| <b>Target</b>          | GNAQ                                                                                      |
| <b>Preservatives</b>   | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol                            |
| <b>Form/Appearance</b> | Liquid                                                                                    |
| <b>Storage</b>         | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                             |
| <b>Note</b>            | For research use only                                                                     |
| <b>Isotype</b>         | IgG                                                                                       |
| <b>Clonality</b>       | Polyclonal                                                                                |
| <b>Uniprot ID</b>      | <b>P50148</b>                                                                             |
| <b>Dilution Range</b>  | ELISA                                                                                     |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                           |