



Product Datasheet

CD4 antibody (orb4830)


Description	CD4 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human CD4. Please contact us for the exact immunogen sequence. The peptide is available as orb374775.
Target	CD4
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	51 kDa
Uniprot ID	P01730 , P05540 , P06332
NCBI	000616 , 000607 , 41 , 11
Entrez	920
Dilution Range	IHC-P: 1:500, IF/ICC: 1:200
Expiration Date	12 months from date of receipt.



WB analysis of human breast cancer (lane...)



WB analysis of human breast tumor (lane ...)



IHC-P image of human thyroid carcinoma t...